# General Protocol for MRK-043-01

MSH, Alpha (Human, Rat, Mouse) - Magnetic Bead RIA kit

(range: 5-640 pg/ml) Sample volume required is 40µl/tube Plasma sample extraction is not required



PHOENIX PHARMACEUTICALS, INC.

#### **INTRODUCTION**

This kit is designed to measure a specific peptide and its related peptides by a competitive radioimmunoassay using Phoenix Particles. It is intended for in vitro protocols only. The Anti-Peptide Magnetic Beads have been conjugated with antiserum against a synthetic form of the peptide.

#### CONTENTS:

1.	RIA buffer, 100ml
2.	Standard peptide, 12.8 µg (lyophilized powder)
3.	Anti-peptide Magnetic Beads, 2.6ml (liquid suspension)
4.	BSA Magnetic Beads for negative control, 0.2ml (liquid suspension)
5.	<sup>125</sup> I-peptide, 1.5 µCi (lyophilized powder)
6.	4x Wash Buffer (2 bottles, 50ml each)
7.	Positive Control (lyophilized powder)
8. Note:	Instructions, 1 booklet and 20 µl pipette tips (12 pieces) Use 20 µl pipette tips for Magnetic Bead aliquoting

Note: Phoenix Pharmaceuticals, Inc. guarantees that its products conform to the information contained in this publication. The purchaser must determine the suitability of the product for their particular needs and establish optimum sample concentration.

## **Storage**

This kit contains reagents sufficient for 125 RIA tubes. Refrigerate all reagents (except <sup>125</sup>I-peptide) between 2 to 8°C for short term storage. **DO NOT FREEZE** Phoenix anti-peptide Magnetic beads and BSA Magnetic beads. <sup>125</sup>I-peptide expires in approximately 6 weeks. Store at -20°C upon receipt. We **strongly recommend** that this kit be used as soon as possible upon receiving it. All solutions should be used on the same day as rehydration.

## **GENERAL INFORMATION**

The assay is based upon the competition of <sup>125</sup>I-peptide and unlabeled peptide (standard or unknown) binding to the limited quantity of antibodies specific for peptide in each reaction mixture. As the quantity of standard or unknown sample in the reaction increases, the amount of <sup>125</sup>I-peptide able to bind to the antibody is decreased. By measuring the amount of <sup>125</sup>I-peptide bound as a function of the concentration of peptide (in standard reaction mixtures), it is possible to construct a "standard curve" from which the concentration of peptide in the unknown sample can be determined. The assay requires one overnight incubation, so plan accordingly.

## Assay Conditions

Plasma, serum, culture media, tissue homogenate, CSF, urine or any biological fluid can be assayed as long as the level of peptide in the sample is high enough for the sensitivity of the kit to detect.

Blood Collection: See page 10.

Tissue Extraction Method: Visit **www.phoenixpeptide.com** and click on the link, "Sample Preparation," for more information.

#### GENERAL PROCEDURE FOR UTILIZATION OF THE RIA KIT:

- 1. Dilute each bottle of 4x Wash Buffer with 150ml of distilled water. This buffer will be used to wash away non-specific binding for the magnetic beads complex after its incubation with sample and tracer.
- 2. Reconstitute the standard peptide with 1ml of RIA buffer. Mix well and store on ice.
- Note: Before adding buffer, carefully examine the eppendorf tube containing the standard. During shipping, part or all of the lyophilized standard may have come loose from the bottom of the tube causing it to stick to the cap or walls of the tube. Gently tap or centrifuge the tube to dislodge powder from the cap or walls. Carefully open the tube and add buffer.

After adding the RIA buffer, vortex for approximately 10 minutes until ALL the peptide powder is completely dissolved. For hydrophobic and hard-to-dissolve peptides, longer vortexing may be required.

- 3. Place the anti-peptide Magnetic Bead solution and BSA Control Magnetic Bead solution on a plate shaker for 20 minutes to resuspend (350rpm, room temperature). If homogeneous suspension is not acheived after 20 minutes, place both bottles in an ultrasonic sonicator, then place on shaker for 20 additional minutes. A homogeneous suspension will not exhibit any bead aggregation at the bottom of the bottle. While the bead solution is still on the shaker (at 350rpm), pipette 20 µl of homogeneous solution to the appropriate tubes (see Table 2).
- 4. Reconstitute the Positive Control with RIA buffer (check QC Data Sheet). Mix well and store on ice.
- 5. Reconstitute samples with RIA buffer (we cannot ensure success with other buffers since they have not been tested).

Tube	Sample	<b>RIA Buffer</b>	Std. Conc.
Stock	Powder	1.0ml	12.8 μg/ml
0	10 µl of Stock	990 μl	128,000 pg/ml
А	5 μl of 0	995 μl	640 pg/ml
В	500 µl of A	500 μl	320 pg/ml
С	500 µl of B	500 μl	160 pg/ml
D	500 µl of C	500 μl	80 pg/ml
Е	500 µl of D	500 μl	40 pg/ml
F	500 µl of E	500 μl	20 pg/ml
G	500 µl of F	500 μl	10 pg/ml
Н	500 µl of G	500 µl	5 pg/ml

6. Prepare dilutions of the standard as in Table 1:

Table 1: Standard Dilutions

- 7. Set up initial RIA reactions (see Table 2) in 12 x 75 mm polystyrene tubes.
  - a) Number tubes TC-1, TC-2, NSB-1, NSB-2, TB-1, TB-2 and #7 #22 for the standards.
  - b) Number tubes #23, #24 for the positive control.
  - c) Number tubes #25 up to #125 for the unkown samples.
  - d) Pipette 140 µl of RIA buffer and 20 µl of BSA Control Magnetic beads into each NSB tube.
  - e) Pipette 140 μl of RIA buffer and 20 μl of anti peptide Magnetic beads into each TB tube.
  - f) Pipette 100 μl of RIA buffer,20 μl anti-peptide Magnetic beads and 40 μl of standards H through A into duplicate tubes #7-#22.
- Note: Reverse the order of preparation so that the concentration increases as the number of the tube increases. For example: Pipette 40 µl of standard H into tubes #7 & #8.
  - g) Pipette 100 µl RIA buffer, 20 µl anti-peptide Magnetic beads and 40 µl of positive control in tubes #23 & #24.
  - h) Pipette 100 μl RIA buffer, 20 μl anti-peptide Magnetic beads and 40 μl of unknown sample into duplicate tubes: tube #25 and up.

- i) Vortex the contents of each tube.
- j) Cover and incubate all tubes for 4 hours at room temperature on an orbital shaker at 350rpm.
- Reconstitute the <sup>125</sup>I-peptide with 12.5ml of RIA buffer and mix well to make tracer solution. Check the concentration of this tracer solution and adjust it with RIA buffer until the concentration (total activity) is 10,000-16,000 cpm/100 μl.
- 9. Add 100  $\mu$ l of the tracer solution to each tube.
- 10. Vortex the contents in each tube.
- 11. Cover and incubate all tubes for 16-24 hours at 4°C.

Tube	Contents	RIA Buffer (ADD FIRST)	Magnetic Beads (ADD SECOND)	Std or Samples (ADD THIRD)	<sup>125</sup> I Peptide (ADD FOURTH)
TC-1 &2	Total Counts				100 µl
NSB-1 & 2	Non-specific binding	140 µl	20 µl BSA		100 µl
TB-1 & 2	Total binding	140 µl	20 µl Anti- peptide		100 µl
7, 8	H Standard	100 µl	20 µl A-P	40 µl	100 µl
9, 10	G Standard	100 µl	20 µl A-P	40 µl	100 µl
11, 12	F Standard	100 µl	20 µl A-P	40 µl	100 µl
13, 14	E Standard	100 µl	20 µl A-P	40 µl	100 µl
15, 16	D Standard	100 µl	20 µl A-P	40 µl	100 µl
17, 18	C Standard	100 µl	20 µl A-P	40 µl	100 µl
19, 20	B Standard	100 µl	20 µl A-P	40 µl	100 µl
21, 22	A Standard	100 µl	20 µl A-P	40 µl	100 µl
23, 24	Positive Control	100 µl	20 µl A-P	40 µl P.C	100 µl
25, 26	Sample 1	100 µl	20 µl A-P	40 µl	100 µl
27, 28	Sample 2	100 µl	20 µl A-P	40 µl	100 µl
Etc.	Etc.	100 µl	20 µl A-P	40 µl	100 µl

Table 2: Contents before Overnight Incubation

- 12. Shake all tubes for one minute at 350 rpm before adding Wash Buffer.
- 13. Add 2000  $\mu$ l of Wash Buffer to each tube except TC tubes. Make sure that beads are suspended in solution; i.e., no pellet is observed at the bottom of any tube. Incubate all tubes at room temperature for 20 minutes. Centrifuge all tubes (except the TC tubes) at 3,000 rpm (approx. 1700 x g) for 20 minutes at 4°C.
- 14. Carefully aspirate **ALL** the supernatant (without touching the pellet) **immediately** following centrifugation (do not decant as the pellet might be lost or excess liquid could be left). **DO NOT ASPIRATE THE TC TUBES.**
- 15. Use a  $\gamma$ -counter to count the cpm of the pellets.
- 16. If results obtained have <15% CV in duplicate tubes, proceed with calculation of results, and skip Steps 17-18.
- Note: If results show >15% variance, this indicates presence of high non-specific binding between samples and beads. In this case, proceed with the additional wash procedure detailed in Steps 17-18.
- If additional wash is necessary, add 1700 μl of Wash Buffer to each tube except TC tubes and vortex the contents in each tube. Incubate all tubes at room temperature for 20 minutes. Centrifuge all tubes (except the TC tubes) at 3,000 rpm (approx. 1700 x g) for 20 minutes at 4°C.
- 18. Repeat Steps 14-15.

## **CALCULATIONS:**

- 1. Using cpm, calculate the average **NSB** and label this as **NSB**.
- 2. Using cpm, calculate the average **TB** and label this as **TB**.
- 3. To find  $B_0$  use the following equation:  $B_0 = TB-NSB$
- 4. To determine the  $B/B_0$  (%) for paired standards and unknown samples use the following calculation:
  - a) Example for standard H:  $B/B_0$  (%) = (Avg. cpm Std. H) - (NSB) x 100%  $B_0$
  - b) Standards G through A (tubes #9-#22), Positive Control (tubes #23 & #24) and the unknown samples (tubes #25 up to #125) are handled as shown above for standard H.
- 5. Examples of tabulated data:

Tube	Samples	Peptide	Average	<b>B</b> / <b>B</b> <sub>0</sub> (%)
			cpm	
TC-1,2			11,170	
NSB-1,2			273	
TB-1,2		0 pg/ml	1,296	100
7,8	H Standard	5 pg/ml	1,267	97.2
9,10	G Standard	10 pg/ml	1,260	96.4
		▲ 		
21,22	A Standard	640 pg/ml	440	16.3
23,24	Positive Control	?	1,096	80.4
25,26	Sample 1	?	976	68.7
27,28	Sample 2	?	1,218	92.3

Table 3: Tabulated Data After Calculation

Total Count (Total activity) (cpm) = 11,170 cpm NSB = 273 cpm TB = 1,296 cpm  $B_0 = 1,296$  cpm - 273 cpm = 1,023 cpm

- 6. On semi log graph paper, plot  $B/B_0$  (%) (in decimal scale ) versus the standard peptide concentrations (in log scale).
  - a) Label the concentrations of standard H through A (5-640 pg/ml) on the X-axis (log scale).
  - b) Label  $B/B_0$  (%) (0 to 100%) on the Y-axis (decimal scale)
  - c) Plot B/B<sub>0</sub> (%) for each standard concentration directly above its X-axis designation.
  - d) Draw the "Best-Fit" curve.
- 7. Determination of the concentration of peptide in unknown samples.
  - a) Using B/B<sub>0</sub> (%) calculated for each unknown sample, read across the graph to the point of intersection with the "Best-Fit" curve.
  - b) The corresponding X-axis coordinate is equivalent to the concentration of peptide (pg/ml) in the assayed sample.
  - c) To calculate the amount of peptide in the original sample, multiply the concentration of the assayed sample by any dilution factor used to prepare the sample.
- Conversion of units: 100 pg/ml x [1000 ÷ Mol. Wt.] = fmole/ml (PMole/L, pM)

# SUMMARY OF ASSAY PROTOCOL

Add anti-peptide Magnetic beads with RIA buffer and samples / standards

Vortex, then incubate on shaker for 4 hours at room temperature

Add <sup>125</sup>I-peptide and vortex.

Incubate for 16-24 hours at 4°C

Vortex/Shake all tubes and add Wash buffer.

Incubate 20 min at room temperature. Centrifuge 20 min at 1,700 x g

Aspirate off the supernatant (except TC tubes)

Count assay tubes and calculate results.

Continue ONLY IF NEEDED. Add Wash buffer.

Incubate 20 min at room temperature. Vortex and centrifuge 20 min at 1,700 x g

Aspirate off the supernatant (except TC tubes)

Count assay tubes and calculate results.

PAGE 9

#### Blood Withdrawal:

Collect blood samples into Lavender Vacutainer tubes (Cat. No. VT-6450) which contain EDTA. Each tube can collect 7ml of blood/tube. Gently rock the Lavender Vacutainer tubes several times immediately after collection of blood for anti-coagulation. Transfer the blood from the Lavender Vacutainer tubes to centrifuge tubes containing aprotinin (Cat. No. RK-APRO) (0.6 TIU/ ml of blood) and gently rock several times to inhibit the activity of proteinases. Centrifuge the blood at 1600 x g for 15 minutes at 4°C and collect the plasma. Plasma kept at -70°C is stable for up to one month.

#### **REFERENCES:**

- Berson, S.A. and Yalow, R.S. Kinetics of reaction between insulin and insulin binding antibody. J. Clin. Invest 36:873, (1957).
- Patrono, C. and Peskar, B.A., (eds) Radioimmunoassay in basic and clinical pharmacology. Heidelberg, Springer-Verlag, (1987).
- Reuter, A., Vrindts-Gevaerts, Y., Meuleman-Gathy, R., Joris, J., Chretien, M. and Franchimont, P. A Radioimmunoassay for Beta-Endorphins. (BETA-END) and (BETA-LPH) in Plasma. Horm Res 25:236, (1987).
- Dwenger, A. Radioimmunoassay: An Overview, J.Clin. Biochem. 22: 883, (1984)
- Wang, Y.N., Chou J., Chang, D., Chang, J.K., Avila, C. and Romero, R. Endothelin-1 in Human Plasma and Amniotic Fluid. In Endothelin-Derived Contracting Factors, edited by G. Rubanyi and P. Vanchoutte, Karger, Basel, pg. 143, (1990).

**CAUTION**: SOME REAGENTS IN THIS KIT CONTAIN SODIUM AZIDE WHICH MAY REACT WITH LEAD AND COPPER PLUMB-ING TO FORM EXPLOSIVE METAL AZIDES. FLUSH WITH LARGE VOLUMES OF WASTE DURING DISPOSAL.

## INSTRUCTIONS FOR POSSESSION, HANDLING AND USE OF RADIOACTIVE MATERIAL

This radioactive material shall only be received, acquired, possessed and used by physicians and veterinarians in clinical laboratories or hospitals for in vitro laboratory tests. Its use should not involve internal or external administration of the material and radiation therefrom to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations and general license requirements of the U.S. Nuclear Regulatory Commission or of a State with which the Commission has entered into an agreement for the exercise of regulatory authority.

### Precautions in Handling Radioactive Material:

The user should store the by-product material, until used, in the original shipping container or in a container providing equivalent radiation protection. There should be no drinking, eating or smoking while radioactive material is being handled. Hands should be covered with gloves during, and thoroughly washed after the handling of radioactive material. When handling radioactive material do not pipette by mouth. Spills must be quickly and thoroughly cleaned up. Surfaces involves should be washed with an alkali detergent (alconox or the equivalent). Persons under 18 should not be permitted to handle radioactive material or enter radioactive areas.

#### Disposal:

Used radioactive test solutions must be disposed of by flushing down a laboratory sink drain with copious amounts of water. Radioactive waste should be disposed of in compliance with Federal, State, and Local Government regulations.

THIS PACKAGE CONFORMS TO THE CONDITIONS AND LIMITA-TIONS SPECIFIED IN 49 CFR173.421 FOR EXCEPTED RADIOAC-TIVE MATERIAL LIMITED QUANTITY, N.O.S. UN2910.

## **Notes**

## **Notes**

#### FOR RESEARCH ONLY

#### NOT FOR USE IN DIAGNOSTIC

**PROCEDURES** 

USA Phoenix Pharmaceuticals, inc.

330 Beach Rd. Burlingame, California 94010 Tel: 650-558-8898, 800-988-1205 Fax: 640-558-1686 Info@phoenixpeptide.com www.phoenixpeptide.com Europe Phoenix Europe GmbH

Viktoriastrasse 3-5. D-76133 Karlsruhe Germany Tel: +49 (721) 16 11 950 Fax: +49 (721) 16 11 952 Germany@phoenixpeptide.com